

RRD

removable rigid dressing use and care



Designed to protect your residual limb, the RRD allows for progressive shrinkage and makes it easy to inspect the wound.

The RRD is worn at all times, except for brief removal 2-3 times a day to allow for limb inspection, dressing changes, and air to the limb. Re-apply after each removal.

- As your residual limb decreases in size, additional socks can be added.
- Remove daily for hygiene, wound inspections, and dressing changes.
- The RRD can be worn with a shrinker sock if recommended by your physician.
- Wear while sleeping to provide wound protection and aid in the prevention of knee contractures.
- If needed, use belt to keep the RRD suspended.
- If you are not wearing the RRD, and you fall, you are not protected.

ANYTIME THERE IS SEVERE DISCOMFORT,
OR ANY OTHER SITUATION THAT DOES NOT SEEM RIGHT,
CALL ELIZUR AT 1-844-628-8813.

Review/Updated 3/19