

DONNING

pin lock system prosthesis



Always check your skin prior to donning your prosthesis for any redness or irritation. If you cannot see the end of your limb, use a hand-held mirror or ask for assistance.

1. Inspect the liner for holes or tears, and ensure the pin has not become loose. If you find any of these conditions, contact your prosthetist immediately.
2. Don the locking liner.
 - a. Make sure the pin is straight and in the center of your limb.
 - b. Make sure there is no air trapped in the end of liner.
3. If you are donning the prosthesis from a sitting position, hold the prosthesis on both sides. Center your limb in the socket, and pull the prosthesis on until you hear a click.
4. Once you hear the click, attempt to pull the prosthesis off to ensure that it is locked before standing.
5. When you stand, the residual limb will go further into the socket and you will hear more clicking. You should hear one click at a time until you are completely settled into the socket.
 - a. If you immediately go all the way to the bottom of the socket, it is too loose. You will need to add socks.
 - b. If it is extremely difficult to engage the lock, remove socks and/or the liner, and realign the locking pin.

**ANYTIME THERE IS SEVERE DISCOMFORT,
OR ANY OTHER SITUATION THAT DOES NOT SEEM RIGHT,
CALL ELIZUR AT 1-844-628-8813.**

Review/Updated 2/19