

# SHRINKER

## *compression tube type application*



## DONNING

1. If rough edges of skin, staples, or stitches are still present, you may need to place a sheath (nylon or cotton) over your limb/dressing to protect the incision.
2. Apply the first shrinker layer, making sure any scabbing, rough edges, staples, or stitches are covered and won't be pulled on by the shrinker. Place the upper edge of the shrinker a couple of inches above the knee.
3. Slide the ring up to the bottom edge of your residual limb. Make sure there are no voids for swelling to occur. Note: You may need to pull your sock up throughout the day.
4. Apply (reflect) the second shrinker layer over the first layer. Make sure the second layer is lower than the first for maximum compression.

ANYTIME THERE IS SEVERE DISCOMFORT,  
OR ANY OTHER SITUATION THAT DOES NOT SEEM RIGHT,  
CALL ELIZUR AT 1-844-628-8813.

Review/Updated 3/19